

# miRCURY LNA™ microRNA Array

## Power Labeling kit

### Instruction manual

for product # 208030-A, 208031-A, 208032-A

### Literature citations:

Please refer to miRCURY LNA™ microRNA Array when describing a procedure for publication using this product.

### Patents and Trademarks

Exiqon, LNA™, Hy3, Hy5 and miRCURY™ are registered trademarks of Exiqon A/S, Vedbaek, Denmark. Locked Nucleic Acids (LNA™) are covered by patents and patent applications owned by Exiqon A/S.

### Disclaimer

Products are for research use only and not for diagnostic or therapeutic use. The products may be used only for the buyer's internal research purposes and not for commercial use. The buyer may not resell products in their original or any modified form. The purchase of products does not include or carry an implied right or license for the buyer to use such products in the provision of services to third parties and a license must be obtained directly from Exiqon A/S for such use.

The miRCURY LNA™ Array labeling kit contains the fluorescent dyes Oyster® 656 and Oyster® 556. Oyster® is a registered trademark of Denovo Biolabels GmbH.

All other trademarks are the property of their respective owners.



# Table of contents

<b>Product summary</b> .....	4
miRCURY LNA™ microRNA Array, Power labeling kit content. ....	4
Additional required material. ....	4
Product description .....	5
Storage .....	6
Related products .....	7
Protocol overview. ....	8
<b>Protocol</b> .....	9
Protocol & Notes .....	11
<b>Tips and Trouble shooting</b> .....	13
<b>References</b> .....	17



# Product summary

4

POWER LABELING KIT | Instruction Manual

## miRCURY LNA™ microRNA Array, Power labeling kit content

The miRCURY LNA™ microRNA Array Power labeling kit consists of 8 vials as described in Table 1.

**Table 1**

Kit Components	Amount supplied	Tube Label	
Calf Intestine Phosphatase (CIP)	20 µL	White	●
CIP buffer	20 µL	Brown	●
Hy5™ fluorescent label	12 or 24 reactions§	Blue	●
Hy3™ fluorescent label	12 or 24 reactions§	Red	●
Labeling enzyme	48 µL	Yellow	●
DMSO	100 µL	Transparent	○
Labeling buffer	150 µL	Orange	●
Nuclease-free water	500 µL	Green	●

§ Product # 208032-A for dual colour labeling contain one tube with the Hy5™ dye and one tube with the Hy3™ dye. The kits for single colour labeling (product # 208030-A and 208031-A) contains two tubes with the Hy3™ dye. Each tube contains sufficient label for 12 labeling reactions.

## Additional required material

- Total RNA containing the small RNAs.
- RNase-free 0.5 mL microcentrifuge tubes.
- Microcentrifuge.
- Incubator or heating block.
- Nuclease-free water.
- miRCURY LNA™ microRNA Array microarray kit.



## Product description

The miRCURY LNA™ microRNA Array Power labeling kit will label RNA molecules with a single fluorophore per molecule. All small RNAs are uniformly labeled. The miRCURY LNA™ microRNA Array labeling kit will label all miRNAs, including those from animals, plants as well as viruses. The labeled RNA molecules may be hybridized to the miRCURY LNA™ pre-spotted microarrays.

The miRCURY LNA™ microRNA Array Power labeling system follows a simple and fast 2-step protocol. The first step includes a Calf Intestinal Alkaline Phosphatase for removal of 5'-phosphates from terminal of the microRNAs. In the second step, a fluorescent label is attached enzymatically to the 3'-end of the miRNAs in the total RNA sample. This is followed by an enzyme inactivation step after which the sample is ready for hybridization.

### Spike-in miRNA Controls

The miRCURY LNA™ microRNA Array Spike-in kit is included in the miRCURY LNA™ microRNA Array, microarray kits, the miRCURY LNA™ microRNA Array, ready-to-spot probe sets, 208110-A and available as a separate product (product# 208040).

The miRCURY LNA™ microRNA Array Spike-in kit contains 10 different synthetic unlabeled miRNAs in different concentrations. The set can be spiked into an RNA sample prior to labeling and the synthetic spike-in miRNAs will hybridize to corresponding capture probes on the miRCURY LNA™ microRNA Array.

For detailed description please see the instruction manual for the miRCURY LNA™ Array Spike-in miRNA kit.



## Storage

The miRCURY LNA™ microRNA Array Power labeling kit is shipped in plastic foam (polystyrene) containers with dry ice. Though it is recommended to store the kit at  $-20^{\circ}\text{C}$ , exposure to higher temperatures (4 to  $10^{\circ}\text{C}$ ) during shipping does not pose any risk to the enzymes.

The fluorescent labeling dyes are shipped dry for increased stability. Before use, add  $30\ \mu\text{L}$  of RNase-free water, vortex and spin to collect tube content. The fluorescent labeling dyes in the kit should not be subjected to repeated cycles of freeze/thawing. Instead it is recommended to aliquot the dyes and store at  $-20^{\circ}\text{C}$ . For long-term storage of the dyes keep the vial(s) at  $-80^{\circ}\text{C}$ . All vials in the kit should be kept on ice during laboratory work and the dyes should always be kept in the dark protected from light.



## Related products

Exiqon offers a tool kit enabling new discoveries concerning the expression, function, and spatial distribution of miRNAs:

### **miRCURY LNA™ microRNA Array, microarray kit**

Pre-printed miRCURY LNA™ microRNA Array microarray slides, available in pack sizes of 3, 6 and 24 (product # 208000-A, 208001-A, 208002-A, 208100-A, 208101-A, 208102-A).

### **miRCURY LNA™ microRNA Array, ready-to-spot probe set**

Ready-to-spot oligo set for direct printing of arrays, or coupling in bead-based applications (product # 208010-A, 208110-A).

### **miRCURY LNA™ microRNA Array, Spike-in miRNA kit**

Ten different synthetic unlabeled miRNAs in different concentrations. The spike-in miRNA kit will hybridize to corresponding capture probes on the miRCURY LNA™ microRNA Array (product # 208040).

### **miRCURY LNA™ microRNA Array, Hybridization buffer**

5 mL hybridization buffer optimal for miRNA hybridization to the miRCURY LNA™ microRNA Arrays (product # 208022).

### **miRCURY LNA™ microRNA Array, Wash buffer kit**

125 mL salt buffer and 15 mL detergent optimal for wash of miRCURY LNA™ microRNA Arrays. (product # 208021).

### **miRCURY LNA™ microRNA Detection**

For in situ hybridization and northern blotting of all annotated miRNAs.

### **miRCURY LNA™ microRNA Knockdown**

miRNA knockdown probes: determine or confirm miRNA function.

### **miRCURY LNA™ microRNA Real-time PCR**

Quickly and accurately determine miRNA expression using real-time PCR system. Available soon.



## Protocol overview

---

### miRCURY LNA™ microRNA Power Labeling Kit

CIP treatment      Mix: RNA sample  
Spike-In miRNA kit



Labeling reaction      Mix: CIP'ed RNA sample  
Labeling buffer  
Hy3™ or Hy5™  
DMSO  
Enzyme



### miRCURY LNA™ microRNA Array Kit

Mix samples      Mix: Hy3™ labeled sample  
Hy5™ labeled sample  
Hybridization buffer  
Denature sample



Hybridize      Hybridize at 56°C for 16 hours



Stringency wash      Wash 2 min. in buffer A at 56°C  
Wash 2 min. in buffer B at 23°C  
Wash 2 min. in buffer C at 23°C  
Dry slides



Image acquisition      Scan slides (recommended scan at 5µm)  
Download relevant GAL files from  
[www.exiqon.com](http://www.exiqon.com)



# Protocol

See tip  
7

## Preparation of RNA sample

Total RNA should be prepared using a method that preserves small RNA species. When using commercially available kits, please ensure that the total RNA preparation is guaranteed to contain microRNAs. An example of a column purification procedure which preserves small RNAs is the miRNeasy Mini Kit from Qiagen. Procedures that include an acidic phenol chloroform extraction are generally recognized as methods that preserve small RNAs. Subsequent to any RNA extraction procedure using a guanidine isothiocyanate-phenol: chloroform extraction, we recommend a column purification procedure to remove traces of these chemicals that potentially may inhibit labeling. An example of a column purification procedure which preserves small RNAs is the miRNeasy Mini Kit from Qiagen. The purified total RNA should be dissolved in RNase-free water at a concentration of 0,5-1  $\mu\text{g}/\mu\text{L}$ .

It is recommended to assess the integrity of the RNA isolated before proceeding with labeling. This may be performed either on the Agilent Bioanalyzer (RIN values should be above 7) or by denaturing gel electrophoresis. Degraded RNA is not suitable for labeling or for hybridization to microarrays.

The procedure used for labeling of miRNAs in the miRCURY LNA™ microRNA Array labeling kit uses total RNA and enrichment for small RNAs is not necessary. Exiqon has carried out extensive comparisons of slides that were hybridized with labeled miRNA enriched samples versus total RNA samples. We do see some differences comparing miRNA enriched and non-enriched miRNA samples. Due to the risk of losing information we do not recommend miRNA enrichment of the total RNA samples.

The amount of total RNA to be labeled for an array hybridization depends on the miRNA content of the cells or tissue being analyzed as this amount is known to vary. Without prior knowledge of miRNA content in the sample to be analyzed we would recommend to use 0,25-1  $\mu\text{g}$  of total RNA per labeling reaction per slide hybridization. Exiqon has compared miRNA expression pattern using a range of 0,25-5  $\mu\text{g}$  labeled total RNA, obtaining same expression pattern in this range.



## Before starting the experiment

Prior to performing the labeling, dissolve the fluorescent dye(s) by adding 29  $\mu\text{L}$  of nuclease-free water to the tube with the labeling dye, followed by vortexing and a brief centrifugation to collect the content of the tube.

### Important note:

When performing dual color hybridization with Hy3<sup>TM</sup> and Hy5<sup>TM</sup> labeled RNA, it is recommended to use Hy3<sup>TM</sup> and Hy5<sup>TM</sup> labels of the same lot number. Do not mix dyes of different lot numbers.

The miRCURY LNA<sup>TM</sup> microRNA Array, Power labeling kit (product # 208032-A) always contains dyes from the same lot#.

To assure optimal labeling and hybridization, ten spike-in control probes (Spike-in miRNA kit) are supplied in the miRCURY LNA<sup>TM</sup> microRNA Array microarray kit and in the ready-to-spot probe set.

Dissolve the spike-in miRNAs in 30  $\mu\text{L}$  of RNase free water supplied upon receipt. Leave the suspension on ice for 30 minutes to dissolve. Vortex and then spin to collect tube contents. Store the dissolved Spike-in miRNA at  $-20^{\circ}\text{C}$  until use and avoid repeated cycles of freeze/thawing. You may wish to aliquot the dissolved spike-in miRNAs to avoid repeated freeze/thawing. For long-term storage, keep the vial at  $-80^{\circ}\text{C}$ .



# Protocol & Notes

## Protocol

Total handling time: 1 hour

**Step 1**  
Thaw all kit components

Place all kit components on ice and thaw for 15-20 min. Mix thoroughly by vortexing followed by brief centrifugation. Do not thaw or vortex the enzymes. Flick these tubes followed by brief centrifugation.

**Step 2**  
Combine reagents according to Table 2. Mix on ice

Reagents should be combined in an RNase-free microcentrifuge tube and should be mixed by pipetting up and down to ensure that all reagents are mixed thoroughly.

Table 2	Volume (µL)
Total RNA**	2
Spike-in miRNA kit	1
CIP buffer	0.5
CIP enzyme	0.5

\*\* We recommend to use between 0.25 - 1 µg, depending on the miRNA content in the sample. Mastermix preparation. In order to minimize variation between slides it is recommended to prepare master mixes for both the CIP reaction (Table 2) and the labeling reaction (Table 3).

**Step 3**  
Incubate 30 min. at 37° C

Incubate 30 min. at 37° C, using a PCR cycler with heated lid.

**Step 4**  
Incubate 5 min. at 95° C

Stop the enzyme reaction and denature the RNA by incubation at 95° C followed by snap cooling on ice.



**Step 5**

Incubate 2 min. on ice.

Leave on ice for at least 2 min. and up to 15 min. Briefly spin the reaction after incubation on ice.

**Step 6**

Combine reagents listed in Table 3. Mix on ice

Add the reagents listed in Table 3 to the 4  $\mu\text{L}$  CIP reaction from step 5.

Table 3	Volume ( $\mu\text{L}$ )
CIP reaction from step 5	4
Labeling buffer	3
Fluorescent label (Hy3™ or Hy5™)	1.5
DMSO	2
Labeling enzyme	2

Mastermix preparation. In order to minimize variation between slides it is recommended to prepare master mixes for both the CIP reaction (Table 1) and the labeling reaction (Table 2).

**Step 7**

Mix and centrifuge the reagents briefly

Reagents should be mixed by gentle vortexing or by pipetting up and down to ensure that all reagents are mixed thoroughly.

**Step 8**

Incubate at 16° C for 1 hour

Incubate for 1 hour at 16° C, using a PCR cycler with heated lid. Protect the reaction from light.

**Step 9**

Incubate for 15 min. at 65° C

After stopping the labeling procedure, briefly spin the reaction and leave it at 4° C. The labeled sample is now ready for hybridization on the array. Hybridization should preferably occur within 1-2 h.



# Tips and Trouble shooting

**Tip 1** **RNA quality**  
Total RNA should be prepared using a method that preserves small RNA species. When using commercially available kits, you should make sure that the total RNA preparation is guaranteed to contain microRNAs. An example of a column purification procedure which preserves small RNAs is the miRNeasy Mini Kit from Qiagen (see tip 7). Procedures that include an acidic phenol chloroform extraction are generally recognized as methods that preserve small RNAs. Subsequent to any RNA extraction procedure using a guanidine isothiocyanate-phenol:chloroform extraction, we recommend a column purification procedure to remove traces of these chemicals that potentially may inhibit labeling. The miRNeasy Mini Kit from Qiagen can be used for this purpose. It is recommended to assess the integrity of the RNA isolated before proceeding with labeling. This may be performed either on the Agilent Bioanalyzer or by denaturing gel electrophoresis. These methods do not reveal/show miRNAs, but yield an assessment of the overall RNA quality.

**Tip 2** **Enrichment for miRNA**  
The miRCURY™ labeling kit will efficiently label RNA that has been enriched for miRNAs as well as miRNAs in total RNA preparations. It should be noted that miRNAs only constitute a minor fraction (~0.01%) of total RNA and attempts to purify this fraction often results in depletion of the miRNAs or a co-purification of large RNAs. We recommend to label and hybridize the total RNA.

**Tip 3** **Weak hybridization signals**  
Several options to investigate the cause of weak signals have been build into the miRCURY LNA™ microRNA Array and labeling kit. Key components of this investigation involves checking the intensities on the landing lights and on the spike-in miRNAs capture probes. Hy5™ label is sensitive to ozone, thus it is important to work fast during entire procedure. If possible handle labeling and especially scanning of the slide in ozone free environment.

**Tip 4** **Landing lights**  
Landing lights are Hy3™ labeled capture probes spotted directly onto the array. If there are very weak signals from the landing lights this could indicate that the PMT settings of the microarray scanner were too low. Consult the instruction manual for your microarray scanner to determine how to choose the optimal PMT settings and redo the scanning. Alternatively, the slide might



have been placed incorrectly in the scanner. Scanning of the slide in an up-side-down or a backwards orientation will result in very little or no signal. Consult the miRCURY LNA™ microRNA Array manual to identify the side and area of the slide where the array has been spotted. See tip 5 & 6.

**Tip 5****Spike-in capture probe**

The miRCURY LNA™ microRNA Array Power Labeling kit can be used with a spike-in kit, consisting of a mixture of synthetic miRNAs. When spike in to the labeling reaction the synthetic miRNAs will co-label with the sample. Subsequently, the synthetic miRNAs will hybridize to control capture probe on the miRCURY LNA™ microRNA Array. If the signals of all spike-in miRNA capture probes are weak, this may indicate that the labeling reaction was performed under suboptimal conditions e.g. pipetting errors, incubation at too low temperatures, exposure of labeling dyes to excess light, loss of miRNAs during concentration of the labeling reaction, presence of RNases or enzyme inhibiting compounds such as phenol or ethanol residues from the sample preparation. Alternatively, hybridization at too stringent conditions (low salt, high temperature) will also result in low signals. If the signal from the spike-in control capture probe is high but signal from most other capture probes are low, this may indicate low quality of the RNA sample (see above) or that the RNA concentration was lower than anticipated.

**Tip 6****Removal of unincorporated dye**

In contrast to chemical labeling procedures, the labeling dye in the The miRCURY LNA™ microRNA Array Power Labeling kit The miRCURY LNA™ microRNA Array Power Labeling kit will not produce high background on the array slide. Therefore, removal of unincorporated dye is normally not necessary. However if removal of dye is desired, we recommend performing an ethanol precipitation. It is important to work fast due to the Hy5™ sensitivity to ozone.



- Mix the Hy3 and Hy5 labeling reactions before precipitation.
- Add 2.5 $\mu$ L RNase free Sodium acetate (3 M, pH 5.5) to the 25 $\mu$ L labeling reaction + 75 $\mu$ L 99.9% ethanol.
- Incubate the sample at -20° C for 30 minutes.
- Centrifuge 30 minutes at max speed in a cooled centrifuge.
- Remove supernatant and wash with 80% precooled ethanol, by centrifuging 5 min at max speed in a cooled centrifuge.
- Remove supernatant and if necessary speedvac for a few minutes to remove remainder of ethanol.
- Dissolve in a mixture of 4 $\mu$ L DMSO, 15 $\mu$ L H<sub>2</sub>O and 6 $\mu$ L labeling buffer.

## Tip 7

**Sample concentration using RNeasy Mini Kit from Qiagen**

RNeasy Mini Kit from Qiagen (product # 74104). The following protocol has been validated and found to concentrate miRNAs with minimal depletion:

- Add 350  $\mu$ L Buffer RLT to the sample, and disrupt and homogenize immediately (Vortex).
- Add 3.5 volumes of 100% ethanol (1225  $\mu$ L), and mix thoroughly by vortexing. Do not centrifuge. Proceed immediately to step 3.
- Pipet 700  $\mu$ L of the sample, including any precipitate that may have formed, into an RNeasy Mini spin column placed in a 2 mL collection tube. Close the lid gently, and centrifuge for 15 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm). Discard the flow-through.

Repeat third step until the whole sample has been pipetted into the spin column. Discard the flow-through each time.

- Place the RNeasy Mini spin column into a new 2 mL collection tube. Pipet 500  $\mu$ L Buffer RPE into the spin column. Close the lid gently, and centrifuge for 15 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm) to wash the spin column membrane. Discard the flow-through.



**Note**

.....  
Buffer RPE is supplied as a concentrate. Ensure that ethanol is added to Buffer RPE before use (see “Things to do before starting” in the handbook supplied with the RNeasy Mini Kit).

- Pipet another 500  $\mu\text{L}$  Buffer RPE into the RNeasy Mini spin column. Close the lid gently, and centrifuge for 15 s at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm) to wash the spin column membrane. Discard the flow-through and the collection tube.
- Place the RNeasy Mini spin column into a new 2 mL collection tube. Centrifuge at full speed for 1 min.
- Place the RNeasy Mini spin column into a 1.5 mL collection tube. Pipet 25  $\mu\text{L}$  RNase-free water directly onto the spin column membrane. Close the lid gently, and centrifuge for 1 min at  $\geq 8000 \times g$  ( $\geq 10,000$  rpm) to elute the miRNA and total RNA.


If the expected RNA yield is  $>30 \mu\text{g}$ , repeat step 7 with a second volume of RNase-free water. Elute into the same collection tube.



# References

- The microRNA Registry.  
Griffiths-Jones S. NAR, 2004, 32, Database Issue, D109-D111
- miRBase, Wellcome Trust Sanger Institute. <http://microrna.sanger.ac.uk>
- [www.exiqon.com/miRCURY/array](http://www.exiqon.com/miRCURY/array)





**Outside North America**

Exiqon A/S · Bygstubben 9  
DK-2950 Vedbaek · Denmark  
Phone +45 45 660 888  
Fax +45 45 661 888

**North America**

Exiqon Inc. · 14 F Gill Street  
Woburn, MA 01801 · United States  
Phone +1 781 376 4150  
Fax +1 781 376 4152  
Toll free [US] +1 888 miRCURY

[support@exiqon.com](mailto:support@exiqon.com)  
[www.exiqon.com](http://www.exiqon.com)

**EXIQON**  
Seek Find Verify